

QUANTACHROME

Quantachrome
INSTRUMENTS



© 2010 Quantachrome Corporation 0710

RIFFLERS

ROTARY POWDER SAMPLE SPLITTERS

07110-Rev.A

REMOVE PART FROM
THREADS ON THIS SIDE

Physical measurements, such as particle size, surface area, density, etc. often require reduction of large quantities into smaller samples. Therefore it is necessary that the powder sub-sample used for analysis be representative of the initial lot or batch for reproducible and valid data. Quantachrome's rifflers are designed to accurately and reproducibly prepare such representative powder samples from larger laboratory or bulk quantities using the statistically superior method of spin, or rotary, riffling.

ROTARY MICRO RIFFLER™

The Rotary Micro Riffler™ easily and automatically divides samples into eight (8) smaller representative samples for analysis. By repetitive divisions of the collected smaller samples, the desired final sample size can be easily obtained.

The Rotary Micro Riffler™ has two basic components: a circular vibrating bowl and a sample collector. The vibrating bowl is amplitude variable and controls the delivery rate of the initial powder sample. The collector consists of a rotating splitter holding eight test tubes. The size of the test tube may be varied to accommodate the initial sample size. Test tubes afford simple filling, emptying and cleaning.

To accommodate a variety of materials the Rotary Micro Riffler™ has two sets of controls: one for the delivery rate and the second for the rotation rate of the collectors. By varying these two controls, accurate divisions of any powder may be achieved, regardless of particle size or density, in a reasonable time.



To divide a sample, up to 120 cm³ of powder is loaded into the vibrator bowl. The vibrator bowl automatically feeds the collector and the test tubes rotating below the feed chute. After all the material has been delivered into the test tubes, further reduction in sample size may be accomplished by emptying the contents of one or more of the test tubes into the bowl and repeating the procedure until the desired sample size is achieved.

Rotary Micro Riffler Specifications

Vibrator Bowl: Aluminum (Nickle-Plated Available)
Max. Capacity: 120 cm³
Collector Tube Volumes: 4 cm³ and 15 cm³
Max. Particle Size: 2.5 mm
Electrical: 110-120V or 220-240V, 50/60 Hz
Height: 33 cm (13")
Width: 20 cm (7.75")
Depth: 18.5 cm (7.25")
Weight: 7.3 kg (16 lbs.)

SIEVING RIFFLER™

The Sieving Riffler™ is designed to handle larger quantities of powder than the Rotary Micro Riffler™. Additionally, particles larger than a particular size and/or foreign bodies can be excluded from the riffling process by attaching a sieve to the top of the hopper.



The Sieving Riffler™ has three basic components: a cylindrical vibrating hopper, a rotating sample collector table and a replaceable sieve. The collector is divided into eight (8) removable stainless steel compartments. The sieve is mounted on top of the vibrating hopper. The Sieving Riffler™ offers the same controls over collector speed and delivery rate as the Rotary Micro Riffler™.

To divide a sample, the hopper is filled or partly filled with sample. The vibrator amplitude is set to automatically feed the powder into the rotating collectors. After all the material has been delivered, further reduction in sample size may be achieved by emptying the contents of one or more collector back into the hopper for additional splitting. When the amount of sample to be re-riffled is no longer reasonable for the Sieving Riffler™, smaller samples can be obtained by splitting the contents of one compartment on the Rotary Micro Riffler™. The durable stainless steel collectors are appropriate for both laboratory and industrial environments.

Sieving Riffler Specifications

Max. Capacity: 2500 cm³
Collector
Compartment Volume: 300 cm³
Max. Particle Size: 5 mm
Sieve Diameter: 5"
Sieve Material : Brass (some sizes also available in Stainless Steel)
Electrical: 110-120V or 220-240V, 50/60 Hz
Height: 51cm (20")
Width: 51cm (20")
Depth: 51cm (20")
Weight: 29 kg (63 lbs.)

Note: When ordering the Sieving Riffler, please specify the sieve size you require.

MICRONS	2360	2000	1700	1180	1000	850	710	600	500	425	355	300	250	212	180	150	125	106	90	75	63	53	45
MESH	8	10	12	16	18	20	25	30	35	40	45	50	60	70	80	100	120	140	170	200	230	270	325



Quantachrome Instruments' corporate headquarters in Boynton Beach, Florida.

Quantachrome®

Renowned innovator of ideas for today's porous materials community.

For over 40 years, Quantachrome's scientists and engineers have revolutionized measurement techniques and designed instrumentation to enable the accurate, precise, and reliable characterization of powdered and porous materials:

- Adsorption/Desorption Isotherms
- Surface Area Measurement
- Pore Size Distribution
- Chemisorption Studies
- Water Sorption Behavior
- Mercury Porosimetry
- True Solid Density
- Tapped Density

Not only are Quantachrome products the instruments of choice in academia, but the technology conceived and developed by our expert staff is applied in industrial laboratories worldwide, where research and engineering of new and improved porous materials is ongoing. Manufacturers also rely on porous materials characterization technology to more precisely specify bulk materials, to control quality, and to isolate the source of production problems with greater efficiency.

Quantachrome is also recognized as an excellent resource for authoritative analysis of your samples in our fully equipped, state-of-the-art powder characterization laboratory.



Quantachrome Instruments Application Laboratory.

CORPORATE HEADQUARTERS

Quantachrome Instruments

1900 Corporate Drive
 Boynton Beach, FL 33426 USA
 Phone: +1 (800) 989-2476
 +1 (561) 731-4999
 Fax: +1 (561) 732-9888
 E-mail: qc.sales@quantachrome.com
 www.quantachrome.com

CHINA

Quantachrome Representative Office

M806. Jingbao Garden
 183 Andingmenwai Street
 Beijing 100011, China
 Phone: +86 800 810 0515
 +86 10 64400892
 +86 13 801191442
 Fax: +86 10 64400892
 www.quantachrome-china.com

EUROPE

Quantachrome UK Limited

Pale Lane Farm, Pale Lane
 Hartley Wintney
 Hook RG27 8BA, UK
 Phone: +44 (0) 1252819719
 Fax: +44 (0) 1252819901
 www.quantachrome.co.uk

EUROPE

Quantachrome GmbH & Co. KG

Rudolf-Diesel Str. 12
 85235 Odelzhausen, Germany
 Phone: +49 (0) 8134/93240
 Fax: +49 (0) 8134/932425
 www.quantachrome.de



Quantachrome Instruments' quality management system is certified to be in accordance with ISO9001:2000.

WORLDWIDE SALES AND SERVICE

- | | |
|-----------------|----------------------|
| Argentina | Mexico |
| Australia | Middle East |
| Austria | Morocco |
| Bahrain | Netherlands |
| Belarus | New Zealand |
| Belgium | Norway |
| Brazil | Oman |
| Bulgaria | Pakistan |
| Canada | Peru |
| Central America | Philippines |
| Chile | Poland |
| China | Portugal |
| Colombia | Puerto Rico |
| Croatia | Romania |
| Cyprus | Russia |
| Czech Republic | Saudi Arabia |
| Denmark | Singapore |
| Egypt | Slovak Republic |
| Estonia | Slovenia |
| Finland | South Africa |
| France | South Korea |
| Germany | Spain |
| Greece | Sri Lanka |
| Hungary | Sweden |
| India | Switzerland |
| Indonesia | Taiwan |
| Ireland | Tanzania |
| Israel | Thailand |
| Italy | Turkey |
| Jamaica | Ukraine |
| Japan | United Arab Emirates |
| Jordan | United Kingdom |
| Kuwait | Uruguay |
| Latin America | Uzbekistan |
| Latvia | Venezuela |
| Lithuania | Vietnam |
| Malaysia | |

Quantachrome

INSTRUMENTS

Serving Porous
Materials and Powder
Characterization
Needs Since 1968



www.quantachrome.com

Trademarks and registered trademarks are the property of their respective owners.